DIAPATH

Constant Innovation in Anatomic Pathology

INSTRUCTION FOR USE

R001X

EN

IT

FR

ES

DE

Acetone

Pure acetone for analysis

Manufacturer

Diapath S.p.A.

Via Savoldini, 71 - 24057 Martinengo (BG) - Italy

Tel. +39 0363 986411

Fax +39 0363 948000

www.diapath.com

info@diapath.com

EMDN code: W01030799

Code	Format	UDI code
R0013	1 LT	08033895009881

Intended use

Disposable, non-sterile in vitro diagnostic device, it performs the function of a pure reagent for analysis in Anatomical pathology. Device for research and diagnostic purposes. The device, used as a support, is intended for manual use, in conjunction with other IVD manual devices.

The device is used in the operational process of an anatomical pathology lab which, once correctly concluded, allows expert and trained personnel to make a diagnosis on the biological sample sectioned.

The device is reserved for professional use only in medical and veterinary laboratories. The device is not reusable once it has reached its end of life. The device does not come into contact with the patient.

For further information, warnings and precautions for use, please read the safety data sheet.

Description of the reagents and any limitations of use

The device consisting of:

Acetone

CAS No. 67-641

EC No. 200-662-2

Device	Stability with closed product	Stability after opening	Ready to use	Storage temperature
R0013	1 year	1 year	YES	15°C – 30°C

Information concerning interfering substances that may affect the performance of the device

Not applicable

List of materials provided and of specific materials that are necessary but not provided

The device is supplied individually and is ready to use.

Information on the combined use of the device

The device can only be used as is or diluted.

Principle of the Test

The reagent performs the role of solvent or reagent for generic use. The device does not produce clinical results in terms of performance.

Instructions for use*

For use of the reagent, refer to the preparation protocols contained in the literature and internal laboratory protocols.

It is neither self-diagnostic device nor is it intended for decentralised analysis. This device may only be used by qualified and suitably trained personnel, since it is mainly intended for health professionals and operators working in specialised medical facilities. Occupational safety guidelines must be applied in accordance with the applicable regulations. The instruments used for making the diagnosis must be suitable for diagnostic use in the laboratory. A diagnosis must only be made by authorised, qualified and trained personnel.

Interpretation of results

Not applicable

Quality control

Not applicable

Analytical performance characteristics

Not applicable

Conditions for sample collection, handling and preparation

The device is designed for processing histological and cytological samples, previously collected and handled by experienced and qualified personnel.

Storage

The device must be stored at a temperature between 15 and 30 °C, as stated on the label.
Keep the container well closed in a cool, well-ventilated place.

The product is stable when properly stored in unopened packaging. Do not use after the expiry date.
If the product is not stored as instructed, its yield may vary and it will have to be tested by the operator.

Stability when handling the device

If incorrectly handled, the device:

- can spill, releasing toxic components into the work environment, thus exposing the operator to the risk of unintentionally inhaling and being exposed to these substances
- may come into superficial contact with the operator, thus causing skin and/or eye irritation
- can compromise the entire sample processing procedure as the device is underperforming

Instructions for disposal / Safe final disposal of the device

Specific instructions on how to safely dispose of the device after use are found in the local regulations of the country where the device is used.

Dispose of contents in accordance with the local/regional regulations.

For information on the content of the reagents, please refer to the safety data sheet. If the product has expired and/or is unused, dispose of it in accordance with local waste regulations, based on the hazard classification on the label, after checking for any contamination. In some cases, analytical evaluation may be required so as to correctly classify the waste and determine the hazard characteristics to be attributed to it

Warnings, precautions and restrictions when using the device

The device must only be used for its intended use.

For proper handling, it is recommended to:

- Use the device only with gloves and other protective PPE
- Do not swallow
- No smoking: flammable device
- Comply with all safety requirements set out for the work environment in which the device is used

- Maintain a clean work environment

Device malfunction or degradation

Due to its chemical characteristics, the device retains its chemical and physical properties throughout its estimated useful life. The device is not subject to deterioration and does not generate debris due to wear and tear and/or degradation products.

Error	Possible cause	Possible corrective action
Failure to maintain a clean work environment	Contamination of the device and noticeable traces of dirt particles	It is recommended to maintain a clean work environment so as to prevent possible contamination of the sample
	Artefact generation	
Delayed diagnosis or non-diagnosis	Device stored under unsuitable temperature conditions	Symbols and appropriate warnings on the device label
Deliberate use beyond the expiry date or intended use limitation	Non-performing device	Symbols and appropriate warnings on the device label

Warnings from the user, with the need to report any serious incidents that have occurred involving the device

In case of a serious accident, the event involving the device in question should be reported to the manufacturer or to the competent authorities, so that the evaluation of such accidents can be more easily coordinated and appropriate safety corrective actions can be taken.

Legend of symbols used in labelling and in the instructions for use

Pictogram	Pictogram name	Explanatory notes
	Manufacturer	This indicates the manufacturer of the device
	Expiry date	This indicates the date after which the device must not be used
	Batch number	This indicates the manufacturer's batch required for batch identification
	Catalogue code	This indicates the code assigned by the manufacturer for identification
	Distributor	This indicates the company that distributes the medical device
	Non-sterile	This indicates that the device has not undergone a sterilisation process
	Keep away from light	This indicates that the device needs to be protected from light sources.
	Temperature limit	This indicates the temperature limit that the medical device can be safely exposed to

	Do not re-use	This indicates a device that is to be used once only, disposable
	Consult instructions for use	This indicates that the user needs to consult the instructions
	Attention	This indicates that the user needs to be careful when operating the device or that the situation requires greater awareness in order to avoid unwanted consequences
	In vitro diagnostic medical device	This indicates a device that is intended to be used as an in vitro diagnostic medical device
	CE Marking	The CE mark certifies that the product complies with the applicable EU directives
	Sufficient content for <n> tests	This indicates the number of tests that can be performed with the device
	Date of manufacture	This indicates the date of manufacture of the medical device
	Unique Device Identification	This indicates a carrier containing the fixed and variable information related to the UDI

Substance/mixture and hazard identification

Substance or mixture	Hazard pictograms	Indication of the hazard	Cautionary advice
Acetone		H225 Highly flammable liquid and vapour	P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. P233 Keep container tightly closed. P261 Avoid breathing dust/fume/gas/mist/vapours/spray P280 Wear protective gloves/protective clothing/Eye protection/face protection. P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. P312 Call a POISON CENTRE or doctor/physician if you feel unwell. P370+P378 In case of fire: use ... to extinguish... P403+P235 Store in a well ventilated place. Keep cool. EUH066 — «Repeated exposure may cause skin dryness or cracking»
		H319 Causes serious eye irritation H336 May cause drowsiness or dizziness	

Bibliography

- “Laboratory methods in histotechnology. Stain methods of the Arm Forces Institute of Pathology” Prophet, Mills, Arrinton, Sobin. Editor McGRAW Hill. Washington D.C. 1968
- “Tecniche di anatomia patologica” M. Melis, F. Carpino, U. Di Tondo. Edi Ermes. 1989

- “Theory and practice of Histological Techniques”. Edited by Bancroft JD and Gamble M. Churchill Livingstone. New York. 2002

Table of revisions

Rev.	Date	Purpose of the revision
1st issue	03/03/2023	Alignment with the new Regulation 746/17

(*) If the request is due to a change, please provide a detailed description of the reason for revising the document and the changes that have been made since the last revision; if available, please enclose recalls, non-conformities, NC CA PA statistics, etc.

You can find this document in the reserved area of the www.diapath.com website. If you want to download the document and keep up to date on the latest available version, you will need to log in after registering.

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